

GC-MS Sample Preparation Protocol for Organic Acids in Urine

1. Scope

This document applies to the preparation and derivatization of urine samples for organic acid analysis by GC-MS.

2. Objective

Illustrate the protocols used for the extraction and derivatization of urinary organic acids for GC-MS analysis.

3. Materials

- 3.1. Urine
- 3.2. HPLC water
- 3.3. 75 g/L (75 mg/mL) methoxyamine hydrochloride solution in pyridine
- 3.4. BSTFA reagent
- 3.5. 3.64 mM cholesterol as internal standard solution
- 3.6. QC synthetic mixture (organic acid mixture, 5 mM)
- 3.7. Ethyl acetate
- 3.8. Hexane

4. Equipment

- 4.1. 1.5 mL Eppendorf tubes
- 4.2. Centrifuge
- 4.3. Adjustable pipet with tips
- 4.4. GC sample vials with 400 uL inserts
- 4.5. Speed Vac Concentrator
- 4.6. Vial Incubator

5. Extraction and Derivatization Protocol

- 5.1. Preparation of blank: 200 μL HPLC water and 40 μL methoxyamine HCl in 2mL glass vials.
- 5.2. Preparation of QC: 100 μL QC mix (200 μM) and 100 μL HPLC water, and 40 μL methoxyamine HCl in 2 mL glass vials.
- 5.3. Pipet samples (200 μL urine) and 40 μL (75 g/L in H_2O) methoxyamine HCl into 2 mL glass vials.
- 5.4. Incubate samples at 60 $^\circ\text{C}$ for 30 minutes.
- 5.5. Transfer samples to 1.5 mL Eppendorf tubes.
- 5.6. Add 20 μL of internal standard (cholesterol, 3.64 mM), and 600 μL of ethyl acetate. Vortex thoroughly for 1 minute. Spin samples at 10000 RPM for three minutes.
- 5.7. Take 500 μL of the supernatant, and put into a new 2 mL glass vial.
- 5.8. Add 600 μL of ethyl acetate to the Eppendorf tube. Vortex thoroughly for one minute. Spin samples at 10000 RPM for three minutes.
- 5.9. Take 500 μL of the supernatant, and put into the 2 mL glass vial containing the previous supernatant (combine the supernatant from the two extractions).
- 5.10. Evaporate samples to dryness under nitrogen with heat (35 $^\circ\text{C}$)
- 5.11. Add 160 μL of Hexane and 40 μL of BSFTA.
- 5.12. Incubate samples at 70 to 90 $^\circ\text{C}$ for 15 minutes.
- 5.13. Transfer samples to 250 μL insert (you can use the same vial).
- 5.14. Refrigerate samples until analysis.
- 5.15. First analyze Hexane, Alkane Standard solution, QC mixture, and Blank using GC-MS. If they all look good (check peaks and Ribitol signal), then proceed to run samples.
- 5.16. In one sequence, set \sim 20 samples to run. Must run QC and Hexane in every 10 samples. And also run QC and Hexane at the end of each sequence.